

Human PPIH/CYP20/CYPH (His tag) recombinant protein



Catalog Number: 519512

General Information

Gene Name Synonym

Peptidyl-Prolyl Cis-Trans Isomerase H, PPIase H, Rotamase H, Small Nuclear Ribonucleoprotein Particle-Specific Cyclophilin H, CypH, U-snRNP-Associated Cyclophilin SnuCyp-20, USA-CYP, PPIH, CYP20, CYPH

Protein Construction

A DNA sequence encoding the Human PPIH/CYP20/CYPH Met1-Met177 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human PPIH/CYP20/CYPH (His tag) recombinant protein