

# Human IGF2BP2/IMP2/VICKZ2 (His tag) recombinant protein



Catalog Number: 519522

## General Information

### Gene Name Synonym

Insulin-Like Growth Factor 2 mRNA-Binding Protein 2, IGF2 mRNA-Binding Protein 2, IMP-2, Hepatocellular Carcinoma Autoantigen p62, IGF-II mRNA-Binding Protein 2, VICKZ Family Member 2, IGF2BP2, IMP2, VICKZ2

### Protein Construction

A DNA sequence encoding the Human IGF2BP2/IMP2/VICKZ2 Met1-Thr220 is expressed with a T7 tag at the N-terminus, 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 12 months at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  as lyophilized powder. Aliquots of reconstituted

samples are stable at  $< -20^{\circ}\text{C}$  for 3 months.

### Formulation

Lyophilized from a  $0.2\ \mu\text{m}$  filtered solution of 20mM PB, 150mM NaCl, pH7.4.

### Usage Guide


### Storage

Lyophilized protein should be stored at  $< -20^{\circ}\text{C}$ , though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at  $4-7^{\circ}\text{C}$  for 2-7 days.

### Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than  $100\ \mu\text{g}/\text{ml}$ . Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

### SDS-PAGE



Human IGF2BP2/IMP2/VICKZ2 (His tag)  
recombinant protein