

# Human SPINK1/PSTI (His tag) recombinant protein



Catalog Number: 519540

## General Information

### Gene Name Synonym

Pancreatic Secretory Trypsin Inhibitor, Serine Protease Inhibitor Kazal-Type 1, Tumor-Associated Trypsin Inhibitor, TATI, SPINK1, PSTI

### Protein Construction

A DNA sequence encoding the Human SPINK1/PSTI Asp24-Cys79 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

## Stability

Store for up to 6 months at -20°C

## Formulation


Supplied as a 0.2 µm filtered solution of 20mM MES, 150mM NaCl, 2mM CaCl<sub>2</sub>, 1mM DTT, 0.05% Brij35, 10% Glycerol, pH 6.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human SPINK1/PSTI (His tag) recombinant protein