

Human UBE2C/UBCH10 (His tag) recombinant protein



Catalog Number: 519550

General Information

Gene Name Synonym

Ubiquitin-Conjugating Enzyme E2 C, UbcH10, Ubiquitin Carrier Protein C, Ubiquitin-Protein Ligase C, UBE2C, UBCH10

Protein Construction

A DNA sequence encoding the Human UBE2C/UBCH10 Met1-Pro179 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human UBE2C/UBCH10 (His tag) recombinant protein