

Human PPP1CC (His tag) recombinant protein



Catalog Number: 519560

General Information

Gene Name Synonym

Serine/Threonine-Protein Phosphatase PP1-Gamma Catalytic Subunit, PP-1G, Protein Phosphatase 1C Catalytic Subunit, PPP1CC

Protein Construction

A DNA sequence encoding the Human PPP1CC Met1-Lys323 is expressed with a 6His tag at the N-terminus, 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 1mM DTT, pH 8.0 .

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human PPP1CC (His tag) recombinant protein