

# Human ATG5/APG5L/ASP recombinant protein



Catalog Number: 519632

## General Information

### Gene Name Synonym

Autophagy protein 5, APG5-like, Apoptosis-specific protein, APG5L, ASP, ATG5,

### Protein Construction

A DNA sequence encoding the Human ATG5/APG5L/ASP Met1-Asp275 is expressed.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

## Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

## Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM Tris, 0.2M NaCl, 1mM DTT, 40% glycerol, pH 8.0 .

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

## SDS-PAGE



Human ATG5/APG5L/ASP recombinant protein