

Human LGALS14/PPL13 (His tag) recombinant protein



Catalog Number: 519642

General Information

Gene Name Synonym

Placental Protein 13-Like, Charcot-Leyden Crystal Protein 2, CLC2, Galectin-14, Gal-14, LGALS14, PPL13

Protein Construction

A DNA sequence encoding the Human LGALS14/PPL13 Met1-Asp139 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris, 100mM NaCl, 1mM DTT, 20% glycerol, pH8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human LGALS14/PPL13 (His tag) recombinant protein