

Human XPNPEP2 (His tag) recombinant protein



Catalog Number: 519645

General Information

Gene Name Synonym

Xaa-Pro Aminopeptidase 2, Aminoacylproline
Aminopeptidase, Membrane-Bound
Aminopeptidase P, Membrane-Bound APP,
Membrane-Bound AmP, mAmP, X-Pro
Aminopeptidase 2, XPNPEP2

Protein Construction

A DNA sequence encoding the Human XPNPEP2
His22-Ala650 is expressed with a 6His tag at the
C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL
method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM
TrisHCl, 150mM NaCl, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human XPNPEP2 (His tag) recombinant protein