

# Human OASL/TRIP14 (His tag) recombinant protein



Catalog Number: 519655

## General Information

### Gene Name Synonym

2'-5'-oligoadenylate synthase-like protein, 2'-5'-OAS-related protein, 2'-5'-OAS-RP, 59 kDa 2'-5'-oligoadenylate synthase-like protein, Thyroid receptor-interacting protein 14, TR-interacting protein 14, TRIP-14, p59OASL, OASL, TRIP14

### Protein Construction

A DNA sequence encoding the Human OASL/TRIP14 Met1-His100 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2  $\mu$ m filtered solution of 20mM PB, 150mM NaCl, pH7.4.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human OASL/TRIP14 (His tag) recombinant protein