

Human OASL/TRIP14 (His tag) recombinant protein



Catalog Number: 519655

General Information

Gene Name Synonym

2'-5'-oligoadenylate synthase-like protein, 2'-5'-OAS-related protein, 2'-5'-OAS-RP, 59 kDa 2'-5'-oligoadenylate synthase-like protein, Thyroid receptor-interacting protein 14, TR-interacting protein 14, TRIP-14, p59OASL, OASL, TRIP14

Protein Construction

A DNA sequence encoding the Human OASL/TRIP14 Met1-His100 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human OASL/TRIP14 (His tag) recombinant protein