

# Human FH recombinant protein



Catalog Number: 519709

## General Information

### Gene Name Synonym:

Fumarate Hydratase Mitochondrial, Fumarase, FH

### Protein Construction:

A DNA sequence encoding the Human FH Ala45-Lys510 is expressed.

**Source:** Human

**Expression Host:** E.coli

## QC Testing

### Activity:

### Purity:

>95% as determined by SDS-PAGE.

### Endotoxin:

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability:

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Predicted N terminal:

## Molecular Mass:

### Formulation:

Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM Tris, pH 8.0.


## Usage Guide

### Storage:

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### Reconstitution:

## SDS-PAGE:



Human FH recombinant protein