

General Information

Gene Name Synonym:

Growth/differentiation factor 11,GDF-11,Bone morphogenetic protein 11,BMP-11

Protein Construction:

A DNA sequence encoding the Human GDF11/BMP11 Asn299-Ser407 is expressed.

Source:Human

Expression Host:Human Cells

QC Testing

Activity:

Purity:

>95% as determined by SDS-PAGE.

Endotoxin:

<1.0 EU per μg protein as determined by the LAL method.

Stability:

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted samples are stable at $< -20^{\circ}\text{C}$ for 3 months.

Predicted N terminal:

Molecular Mass:

Formulation:

Lyophilized from a $0.2\ \mu\text{m}$ filtered solution of PBS,pH7.4.

Usage Guide


Storage:

Lyophilized protein should be stored at $<-20^{\circ}\text{C}$, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at $4-7^{\circ}\text{C}$ for 2-7 days.

Reconstitution:

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than $100\ \mu\text{g}/\text{ml}$. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE:



Human GDF11/BMP11 recombinant protein