

Human CX3CL1/FKN/NTT/SCYD1/A-152E5.2 (His tag) recombinant protein



Catalog Number: 519792

General Information

Gene Name Synonym:

Fractalkine, C-X3-C Motif Chemokine 1, CX3C Membrane-Anchored Chemokine, Neurotactin, Small-Inducible Cytokine D1, CX3CL1, FKN, NTT, SCYD1

Protein Construction:

A DNA sequence encoding the Human CX3CL1/FKN/NTT/SCYD1/A-152E5.2 Gln25-Arg339 is expressed with a 6His tag at the C-terminus.

Source: Human

Expression Host: Human Cells

QC Testing

Activity:

Purity:

>95% as determined by SDS-PAGE.

Endotoxin:

<1.0 EU per µg protein as determined by the LAL method.

Stability:

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Predicted N terminal:

Molecular Mass:

Formulation:

Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Usage Guide


Storage:

Lyophilized protein should be stored at <-20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

Reconstitution:

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE:



Human CX3CL1/FKN/NTT/SCYD1/A-152E5.2 (His tag) recombinant protein