

Human FEN1/RAD2 recombinant protein



Catalog Number: 519866

General Information

Gene Name Synonym:

Flap Endonuclease 1, FEN-1, DNase IV, Flap Structure-Specific Endonuclease 1, Maturation Factor 1, MF1, hFEN-1, FEN1, RAD2

Protein Construction:

A DNA sequence encoding the Human FEN1/RAD2 Met1-Lys380 is expressed.

Source:Human

Expression Host:E.coli

QC Testing

Activity:

Purity:

>95% as determined by SDS-PAGE.

Endotoxin:

<1.0 EU per μ g protein as determined by the LAL method.

Stability:

Store for up to 6 months at -20°C

Predicted N terminal:

Molecular Mass:

Formulation:

Supplied as a 0.2 μ m filtered solution of 50 mM Tris, 50mM NaCl, 1mM DDT, 10% Glycerol, pH 8.0.


Usage Guide

Storage:

Store at <-20°C, minimize freeze-thaw cycles.

Reconstitution:

SDS-PAGE:



Human FEN1/RAD2 recombinant protein