

Human TIGIT/VSIG9/VSTM3 (His tag) recombinant protein



Catalog Number: 519876

General Information

Gene Name Synonym:

T-cell immunoreceptor with Ig and ITIM domains,,VSIG9, VSTM3,TIGIT,V-set and transmembrane domain-containing protein 3,V-set and immunoglobulin domain-containing protein 9

Protein Construction:

A DNA sequence encoding the Human TIGIT/VSIG9/VSTM3 Met22-Pro141 is expressed with a 6His tag at the C-terminus.

Source:Human

Expression Host:Human Cells

QC Testing

Activity:

Purity:

>95% as determined by SDS-PAGE.

Endotoxin:

<1.0 EU per μg protein as determined by the LAL method.

Stability:

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted samples are stable at $< -20^{\circ}\text{C}$ for 3 months.

Predicted N terminal:

Molecular Mass:

Formulation:

Lyophilized from a $0.2\ \mu\text{m}$ filtered solution of PBS,pH7.4.

Usage Guide


Storage:

Lyophilized protein should be stored at $<-20^{\circ}\text{C}$, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at $4-7^{\circ}\text{C}$ for 2-7 days.

Reconstitution:

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than $100\mu\text{g/ml}$. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE:



Human TIGIT/VSIG9/VSTM3 (His tag)
recombinant protein