

Human BTN2A1/BT2.1/BTF1 (His tag) recombinant protein



Catalog Number: 519882

General Information

Gene Name Synonym

Butyrophilin subfamily 2 member A1; BTN2A1;
BT2.1; BTF1

Protein Construction

A DNA sequence encoding the Human
BTN2A1/BT2.1/BTF1 Gln29-Ala248 is expressed
with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL
method.

Stability

Store for up to 12 months at -20°C to -80°C as
lyophilized powder. Aliquots of reconstituted
samples are stable at $< -20^{\circ}\text{C}$ for 3 months.

Formulation

Lyophilized from a $0.2\ \mu\text{m}$ filtered solution of PBS,
 $\text{pH}7.4$.

Usage Guide


Storage

Lyophilized protein should be stored at $<-20^{\circ}\text{C}$,
though stable at room temperature for 3 weeks.
Reconstituted protein solution can be stored at
 $4-7^{\circ}\text{C}$ for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not
mix by vortex or pipetting. It is not recommended
to reconstitute to a concentration less than
 $100\mu\text{g/ml}$. Dissolve the lyophilized protein in
 ddH_2O . Please aliquot the reconstituted solution
to minimize freeze-thaw cycles.

SDS-PAGE



Human BTN2A1/BT2.1/BTF1 (His tag)
recombinant protein