

# Mouse Efnb1/Epl2/Eplg2/Lerk2/Stra1 (Fc & His tag) recombinant protein



Catalog Number: 519904

## General Information

### Gene Name Synonym

Ephrin-B1, EFL-3, ELK ligand, ELK-L, EPH-related receptor tyrosine kinase ligand 2, LERK-2, EFNB1, EFL3, EPLG2, LERK2

### Protein Construction

A DNA sequence encoding the Mouse Efnb1/Epl2/Eplg2/Lerk2/Stra1 Lys30-Ser229 is expressed with a Fc, 6His tag at the C-terminus.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

### Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

## Formulation

Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.

## Usage Guide


### Storage

Lyophilized protein should be stored at <-20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

### Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## SDS-PAGE



Mouse Efnb1/Epl2/Eplg2/Lerk2/Stra1 (Fc & His tag) recombinant protein