

Human MCL1/BCL2L3 (His tag) recombinant protein



Catalog Number: 519909

General Information

Gene Name Synonym

Induced myeloid leukemia cell differentiation protein Mcl-1; Bcl-2-like protein 3; Bcl2-L-3; Bcl-2-related protein EAT/mcl1; mcl1/EAT; MCL1; BCL2L3

Protein Construction

A DNA sequence encoding the Human MCL1/BCL2L3 Met1-Gly327 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM Tris, 200mM NaCl, 2mM DTT, 20% glycerol, pH8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human MCL1/BCL2L3 (His tag) recombinant protein