

Human TNFRSF9/CD137/ILA (Fc tag) recombinant protein



Catalog Number: 519927

General Information

Gene Name Synonym:

CD137, ILA, TNFRSF9, 4-1BB ligand receptor, CDw137, T-cell antigen 4-1BB homolog, T-cell antigen ILA

Protein Construction:

A DNA sequence encoding the Human TNFRSF9/CD137/ILA Leu24-Gln186 is expressed with a Fc tag at the C-terminus.

Source:Human

Expression Host:Human Cells

QC Testing

Activity:

Purity:

>95% as determined by SDS-PAGE.

Endotoxin:

<1.0 EU per μg protein as determined by the LAL method.

Stability:

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted samples are stable at $< -20^{\circ}\text{C}$ for 3 months.

Predicted N terminal:

Molecular Mass:

Formulation:

Lyophilized from a $0.2\ \mu\text{m}$ filtered solution of PBS, pH7.4.

Usage Guide


Storage:

Lyophilized protein should be stored at $<-20^{\circ}\text{C}$, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at $4-7^{\circ}\text{C}$ for 2-7 days.

Reconstitution:

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than $100\ \mu\text{g}/\text{ml}$. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE:



Human TNFRSF9/CD137/ILA (Fc tag)
recombinant protein