

Human SULT4A1/SULTX3 recombinant protein



Catalog Number: 519931

General Information

Gene Name Synonym

Sulfotransferase 4A1, ST4A1, Brain
Sulfotransferase-Like Protein, hBR-STL, hBR-STL-1, Nervous System Sulfotransferase, NST, SULT4A1, SULTX3

Protein Construction

A DNA sequence encoding the Human SULT4A1/SULTX3 Met1-Leu284 is expressed.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human SULT4A1/SULTX3 recombinant protein