

DLX4 antibody

Catalog Number: 110000

Product name

DLX4 antibody

IHC: 1:20-1:200

IF: 1:20-1:200

Specificity

human, mouse, rat; other species not tested.

Antibody description

DLX4 Rabbit Polyclonal antibody. Positive IF detected in HeLa. Positive IHC detected in human lymphoma. Positive WB detected in human heart tissue. Observed molecular weight by Western-blot: 38kd

Preparation

This antibody was obtained by immunization of DLX4 recombinant protein (Accession Number: NM_138281). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:110000(DLX4 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:110000(DLX4 antibody) at dilution of 1:50 (under 10x lens)



Immunofluorescent analysis of HeLa cells, using DLX4 antibody Catalog No:110000 at 1:50 dilution

DLX4 antibody



Catalog Number: 110000

and FITC-labeled donkey anti-rabbit IgG (green). Blue pseudocolor = DAPI (fluorescent DNA dye).