

Anti-GFRA2 antibody



Catalog Number: 106000

Product name

Anti-GFRA2 antibody

Immunogen

[Human GFRA2 \(His Tag\) recombinant protein](#)

Specificity

Human GFRA2

Antibody description

Rabbit monoclonal to GFRA2

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human GFRA2 (rh GFRA2; NP_001486.4; Met1-Ser441).

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid

contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 20-50 µg/mL

Validations



Items	Lanes	A
Sample (whole cell lysate)		HeLa
Sample Volume (µg/lane)		60
Gel		15% SDS-PAGE reducing gel
Recommended Concentration		20-50 µg/ml
Secondary Antibody		DyLight 600 labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.
Developed using Odyssey imaging system.		
Explanation		Predicted band size: 51 kDa Observed band size: 61 kDa Additional bands at: 36 kDa (We are unsure as to the identity of this extra band. It may be a non-specific binding)

GFRA2 Antibody, Rabbit MAb, Western blot