

# Imaprelimab reference antibody



Catalog Number: 177000

## Product name

Imaprelimab reference antibody

## Specificity

Human MUC18 / MCAM / CD146

## Antibody description

Anti-MUC18 / MCAM / CD146 Antibody  
(Imaprelimab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

1×PBS, pH 6.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG

## Applications


ELISA

## Validations



## SDS-PAGE

Anti-MUC18 / MCAM / CD146 Antibody (imaprelimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



## SEC-HPLC

The purity of Anti-MUC18 / MCAM / CD146 Antibody (imaprelimab) is more than 95%, determined by SEC-HPLC.