

ZNF509 antibody



Catalog Number: 117000

Product name

ZNF509 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZNF509 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 85 kDa

Preparation

This antibody was obtained by immunization of ZNF509 recombinant protein (Accession Number: NM_001330625). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:117000(ZNF509 Antibody) at dilution of 1:600