

Anti-Methionine Aminopeptidase 2/p67 Rabbit antibody



Catalog Number: 168000

Product name

Rabbit IgG

Anti-Methionine Aminopeptidase 2/p67 Rabbit antibody

Applications

WB, IHC-P, ICC/IF, FC

Specificity

Human, Mouse, Rat

Dilutions

WB: 1/1000

IHC: 1/50

ICC/IF: 1/50

FC: 1/20-1/100

Antibody description

Methionine Aminopeptidase 2/p67 Rabbit monoclonal antibody

Validations

Preparation

Antigen: A synthetic peptide of human Methionine Aminopeptidase 2/p67

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of Methionine Aminopeptidase 2/p67 in A549,HL-60,U251,U87-MG cell lysates using Methionine Aminopeptidase 2/p67 Rabbit mAb(1:1000 diluted).Predicted band size:53kDa.Observed band size:67kDa.