

# MYO1E Polyclonal Antibody



Catalog Number: 161400

## Product name

MYO1E Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to MYO1E

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 809-1108 of human MYO1E (NP\_004989.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:1000 - 1:2000

## Validations



Western blot - MYO1E Polyclonal Antibody

Western blot analysis of extracts of HeLa cells, using MYO1E antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.