

# NTHL1 antibody

Catalog Number: 113400

## Product name

NTHL1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NTHL1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 34 kDa

## Preparation

This antibody was obtained by immunization of NTHL1 recombinant protein (Accession Number: NM\_002528). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113400(NTHL1 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:113400(NTHL1 antibody) at dilution of 1:50 (under 10x lens)