

# Human YES1/c-Yes (His & GST Tag) recombinant protein



Catalog Number: 501800

## General Information

### Gene Name Synonym

Proto-oncogene c-Yes; p61-Yes

### Protein Construction

A DNA sequence encoding the human YES1 (NP\_005424.1) (Gly 2-Leu 543) was fused with the N-terminal polyhistidine-tagged GST tag at the N-terminus.

### Organism

Human

### Expression Host

Baculovirus-Insect Cells

## QC Testing

### Activity

The specific activity was determined to be 35 nmol/min/mg using Poly(Glu,Tyr) 4:1 as substrate.

### Purity

> 80 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Met

### Molecular Mass

The recombinant human YES1/GST chimera consists of 779 amino acids and has a calculated molecular mass of 88.5 kDa. It migrates as an approximately 75 kDa band in SDS-PAGE under reducing conditions.

### Formulation

Supplied as sterile 20mM Tris, 500mM NaCl, 10% gly, 0.5mM TCEP, pH 8.0

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human c-Yes / YES1 Protein (His & GST Tag)  
SDS-PAGE