

PHAX antibody



Catalog Number: 113800

Product name

PHAX antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PHAX Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of PHAX recombinant protein (Accession Number: NM_032177). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113800(PHAX antibody) at dilution of 1:800