

# Anti-CASP4 antibody



Catalog Number: 175800

## Product name

Anti-CASP4 antibody

## Specificity

Human

## Antibody description

Rabbit polyclonal antibody to CASP4

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human caspase 4/5. aa range:221-270

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB

## Dilutions

WB: 1:500-1:2000

## Validations




Fig1:; Western blot analysis of Caspase 4/5 (p20, Cleaved-Asp270/Asp311) Antibody. The lane on the right is blocked with the Caspase 4/5 (p20, Cleaved-Asp270/Asp311) peptide.