

# Anti-BIN1 Rabbit antibody



Catalog Number: 166800

## Product name

Anti-BIN1 Rabbit antibody

WB: 1/1000

IHC: 1/20

## Specificity

Human, Mouse, Rat

ICC/IF: 1/20

FC: 1/20

## Antibody description

BIN1 Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human BIN1

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications


WB, IHC-P, ICC/IF, FC

## Dilutions

## Validations



Immunohistochemistry of BIN1 in paraffin-embedded Human tonsil using BIN1 Rabbit pAb at



Western blot detection of BIN1 in Rat Brain lysates using BIN1 Rabbit pAb(1:1000 diluted). Predicted band size:65kDa.Observed band size:45-80kDa.