

# ABCB9 antibody



Catalog Number: 107800

## Product name

ABCB9 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ABCB9 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in mouse thymus tissue, HL-60 cells, human brain tissue, human testis tissue, K-562 cells, mouse brain tissue. Observed molecular weight by Western-blot: 72 kDa

## Preparation

This antibody was obtained by immunization of ABCB9 recombinant protein (Accession Number: BC017348). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:107800(ABCB9 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using ABCB9 antibody Catalog No:107800 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).Blue pseudocolor = DAPI (fluorescent DNA dye).