

# C9orf46 antibody



Catalog Number: 108800

## Product name

C9orf46 antibody

## Specificity

Human; other species not tested.

## Antibody description

C9orf46 Rabbit Polyclonal antibody. Positive WB detected in A431 cells, human stomach tissue. Observed molecular weight by Western-blot: 17-20 kd

## Preparation

This antibody was obtained by immunization of C9orf46 recombinant protein (Accession Number: XM\_011517960). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:108800(C9orf46 Antibody) at dilution of 1:1000