

GMIP antibody



Catalog Number: 111100

Product name

GMIP antibody

Immunogen

[Human GMIP Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

GMIP Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 130 kDa

Preparation

This antibody was obtained by immunization of GMIP recombinant protein (Accession Number: NM_016573). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111100(GMIP Antibody) at dilution of 1:600



Immunofluorescent analysis of HepG2 cells using Catalog No:111100(GMIP Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)