

ANO10 antibody

Catalog Number: 108100

Product name

ANO10 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ANO10 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 76 kDa

Preparation

This antibody was obtained by immunization of ANO10 recombinant protein (Accession Number: BC038855). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108100(ANO10 Antibody) at dilution of 1:1000



IP Result of anti-ANO10 (IP:Catalog No:108100, 3ug; Detection:Catalog No:108100 1:500) with mouse brain tissue lysate 3600ug.