

# Anti-GAPDH Mouse mAb (HRP Conjugate)



Catalog Number: 169100

## Product name

Anti-GAPDH Mouse mAb (HRP Conjugate)

## Specificity

Human, Mouse, Rat

## Antibody description

GAPDH Mouse mAb (HRP Conjugate)

## Preparation

Antigen: A synthetic peptide corresponding to residues near the amino terminus of human GAPDH protein.

## Formulation

In PBS(pH 7.4) containing with 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

WB, ELISA

## Dilutions

WB: 1/500—1/5000

## Validations



Western blot detection of GAPDH in Rat brain, CHO-K1 and HUVEC cell lysates using GAPDH Mouse mAb (HRP Conjugate) (1:5000 diluted). Predicted band size: 36KDa. Observed band size: 36KDa.