

Anti-GAPDH Mouse mAb (HRP Conjugate)



Catalog Number: 169100

Product name

Anti-GAPDH Mouse mAb (HRP Conjugate)

Specificity

Human, Mouse, Rat

Antibody description

GAPDH Mouse mAb (HRP Conjugate)

Preparation

Antigen: A synthetic peptide corresponding to residues near the amino terminus of human GAPDH protein.

Formulation

In PBS(pH 7.4) containing with 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB, ELISA

Dilutions

WB: 1/500—1/5000

Validations



Western blot detection of GAPDH in Rat brain, CHO-K1 and HUVEC cell lysates using GAPDH Mouse mAb (HRP Conjugate) (1:5000 diluted). Predicted band size: 36KDa. Observed band size: 36KDa.