

Human Atlastin-3/ATL3 (His Tag) recombinant protein



Catalog Number: 501900

General Information

Protein Construction

A DNA sequence encoding the human ATL3 (Q6DD88) (Met1-Ala445) was expressed with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> (7.5+29.7+61.2) % as determined by SDS-PAGE.

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human ATL3 consists 463 amino

acids and predicts a molecular mass of 52.6 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human ATL3 / Atlastin GTPase 3 Protein (His Tag)
SDS-PAGE