

Anti-TMEM30B antibody



Catalog Number: 175900

Product name

Polyclonal

Anti-TMEM30B antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse

Applications

Antibody description

WB

Rabbit polyclonal antibody to TMEM30B

Dilutions

Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tmem30b. aa range:91-140

Validations

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality




Fig1.; Western Blot analysis of K562 cells using Cdc50B Polyclonal Antibody