

# AF-6 antibody

Catalog Number: 107900

---

## Product name

Polyclonal

AF-6 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

AF-6 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, A431 cells, human brain tissue. Observed molecular weight by Western-blot: 182-210kd

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_005936). Purification method: Antigen affinity purified.

## Validations



## Formulation

0.1M NaHCO<sub>3</sub>, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:107900(MLLT4 antibody) at dilution of 1:500

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality