

CBLC antibody



Catalog Number: 108900

Product name

CBLC antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CBLC Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, A549 cells, COLO 320 cells, HeLa cells, HepG2 cells, L02 cells, rat colon tissue. Observed molecular weight by Western-blot: 52kd

Preparation

This antibody was obtained by immunization of CBLC recombinant protein (Accession Number: NM_012116). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:108900(CBLC antibody) at dilution of 1:500