

Anti-CD8a / Lyt2 antibody



Catalog Number: 101200

Product name

Anti-CD8a / Lyt2 antibody

Immunogen

[Ferret CD8a / Lyt2 \(His Tag\) recombinant protein](#)

Specificity

Ferret CD8a / Lyt2

Antibody description

Rabbit monoclonal to CD8a / Lyt2

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Ferret CD8a / Lyt2 (rF CD8a / Lyt2; ABS50091.1; Met1-Glu186) and conjugated with PerCP under optimum conditions, the unreacted PerCP was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze !Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



CD8a / Lyt2 Antibody (PerCP), Rabbit MAb, Flow cytometric analysis

Analysis of Ferret CD8A expression on ferret splenocytes. Ferret splenocytes were stained with PerCP-conjugated rabbit anti-Ferret CD8A. The histogram were derived from the gated events based on light scattering characteristics of viable cells.