

NKRF antibody

Catalog Number: 113200

Product name

NKRF antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NKRF Rabbit Polyclonal antibody. Positive IP detected in A549 cells. Positive WB detected in Jurkat cells, A431 cells, A549 cells, HepG2 cells, PC-3 cells. Observed molecular weight by Western-blot: 100 kDa

Preparation

This antibody was obtained by immunization of NKRF recombinant protein (Accession Number: NM_001173488). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113200(NKRF antibody) at dilution of 1:1000



IP Result of anti-NKRF (IP:Catalog No:113200, 3ug; Detection:Catalog No:113200 1:1000) with A549 cells lysate 1000ug.