

Visilizumab (Synonyms: MAB HUMANIZED TO CD3+ T CELLS)



Catalog Number: 177200

Product name

IgG2SA

Visilizumab (Synonyms: MAB HUMANIZED TO CD3+ T CELLS)

Applications

ELISA


Specificity

Human CD3

Validations

Antibody description

Anti-CD3 Antibody (Visilizumab)



Preparation

Recombinant expression and purified from CHO cells.

SDS-PAGE


Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Anti-CD3 Antibody (visilizumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.



Clonality

Monoclonal

SEC-HPLC

Ig Type

The purity of Anti-CD3 Antibody (visilizumab) is more than 95%, determined by SEC-HPLC.