

# Visilizumab (Synonyms: MAB HUMANIZED TO CD3+ T CELLS)



Catalog Number: 177200

## Product name

IgG2SA

Visilizumab (Synonyms: MAB HUMANIZED TO CD3+ T CELLS)

## Applications

ELISA


## Specificity

Human CD3

## Validations

## Antibody description

Anti-CD3 Antibody (Visilizumab)



## Preparation

Recombinant expression and purified from CHO cells.

## SDS-PAGE


## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Anti-CD3 Antibody (visilizumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.



## Clonality

Monoclonal

## SEC-HPLC

## Ig Type

The purity of Anti-CD3 Antibody (visilizumab) is more than 95%, determined by SEC-HPLC.