

FAM76A antibody



Catalog Number: 110600

Product name

FAM76A antibody

Immunogen

[Human FAM76A Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FAM76A Rabbit Polyclonal antibody. Positive IHC detected in human placenta tissue, human brain tissue. Positive WB detected in HeLa cells, COLO 320 cells. Observed molecular weight by Western-blot: 38 kDa

Preparation

This antibody was obtained by immunization of FAM76A recombinant protein (Accession Number: NM_152660). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110600(FAM76A antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human placenta using Catalog No:110600(FAM76A antibody) at dilution of 1:50 (under 10x lens)



FAM76A antibody



Catalog Number: 110600

Immunohistochemical of paraffin-embedded human placenta using Catalog

No:110600(FAM76A antibody) at dilution of 1:50 (under 40x lens)