

BCL7B Polyclonal Antibody



Catalog Number: 164300

Product name

BCL7B Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to BCL7B

Preparation

Antigen: A synthetic Peptide of human BCL7B

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:1000

IHC 1:20 - 1:50

Validations



Western blot - BCL7B Polyclonal Antibody

Western blot analysis of extracts of mouse heart, using BCL7B antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.