

Anti-Hsp60 (3G8) Mouse antibody



Catalog Number: 168300

Product name

Anti-Hsp60 (3G8) Mouse antibody

Specificity

Human, Rat, Mouse, Monkey

Antibody description

Hsp60 (3G8) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human HSP60 expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, IHC-P, ICC/IF, FC

Dilutions

WB: 1/500 - 1/2000


IHC: 1/200 - 1/1000

ICC: 1/200 - 1/1000

FC: 1/200 - 1/400

ELISA: 1/10000

Validations



Western blot detection of Hsp60 in Jurkat, HepG2 and HeLa cell lysates using Hsp60 mouse mAb (dilution 1:500). Predicted band size: 61KDa. Observed band size: 61KDa.