

# RILP antibody



Catalog Number: 114700

## Product name

RILP antibody

## Immunogen

[Human RILP Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

RILP Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HepG2 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 41-45 kDa

## Preparation

This antibody was obtained by immunization of RILP recombinant protein (Accession Number: NM\_031430). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114700(RILP Antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells using Catalog No:114700(RILP Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)