

Anti-UGP2 antibody



Catalog Number: 176700

Product name

Anti-UGP2 antibody

Specificity

Human, Mouse, Rat, Zebrafish

Antibody description

Rabbit polyclonal antibody to UGP2

Preparation

This antigen of this antibody was peptide

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:1,000-1:2,000

Validations




Fig1: Western blot analysis on cell lysates using anti-UGP2 rabbit polyclonal antibodies.